

# BT-200 Fetal Doppler

## HI-bebe

**HI-bebe is showing you what fetus is saying.**

HI-bebe is pocket size doppler for Obstetricians, midwives and pregnant women.

- Compact and light
- High-Sensitivity Doppler probe
- Ergonomic design
- Large-Display & clean sound
- Easy to use
- Long time continuous use (6 hours)



Sound version



# BT-200 Fetal Doppler **Hi-bebe**

## Specifications

|                           |                                                                                                       |                         |
|---------------------------|-------------------------------------------------------------------------------------------------------|-------------------------|
| Display                   | <ul style="list-style-type: none"> <li>• LCD Back light</li> <li>• Singnal quality display</li> </ul> | • Low battery display   |
| Heart rate range          | • 50~240 bpm                                                                                          |                         |
| FHR Accuracy              | • $\pm 2\%$                                                                                           |                         |
| Ultrasound Frequency      | • 2 MHz                                                                                               |                         |
| Ultrasound Intensity      | • $<10 \text{ m W/cm}^2$                                                                              |                         |
| Speaker output power      | • 1.2 W (Built in Loudspeaker)                                                                        |                         |
| Audio Output for Earphone |                                                                                                       |                         |
| Auto Shut off Sound       | • 1 Minute                                                                                            |                         |
| Auto Shut off Power       | • 5 Minutes                                                                                           |                         |
| Power                     | • 1.5 V Battery X 2 (Type : AA)                                                                       |                         |
| Sensitivity               | • 10-12 Weeks onward                                                                                  |                         |
| PC Interface              | • Sound card (BCM-220)                                                                                |                         |
| Battery Life              | • 360 Minutes                                                                                         |                         |
| Accessories               | <ul style="list-style-type: none"> <li>• Carrying case</li> <li>• Gel (60ml)</li> </ul>               | • Battery (1.5V AA X 2) |

## Waterproof (Optional)

## PC Software

### BCM220

- Record and play fetal heart sound (.wav format)

