

M3A Vital Signs Monitor

EDAN M3A Vital Signs Monitor provides essential measurement of patient's SpO₂, NIBP and Temperature. It offers an effective solution and brings high performance vital signs monitoring at your finger tips.

12,000

NIBP Measurement
Review

100 h

Trend Review

- 3.5" high resolution color TFT-LCD screen & LED display
- Lightweight & portable design
- User-friendly interface
- Flexible configurations to meet different clinical needs
- Auto/Manual/Continuous/Average BP mode
- Real-time measurements & trend display
- Built-in rechargeable Lithium-ion battery for 17 hours of continuous work
- Compatible with PatientCare Viewer PC data management software
- Nurse call
- Suitable for adult, pediatric and neonatal patients
- Barcode scanner support (optional)
- Ethernet connection for bi-directional communication with MFM-CMS central monitoring system (optional)
- Wi-Fi support with optional AP



Web Viewer



MFM-CMS




Mounting Solutions



Wall Mount



Rolling Stand



Bed Rail Clamp

Specification

Configuration	SpO ₂	NIBP	Covidien Quick Temp
SpO ₂ only	✓	—	—
NIBP only	—	✓	—
SpO ₂ +NIBP	✓	✓	—
NIBP+Covidien Quick Temp	—	✓	✓
SpO ₂ +NIBP+Covidien Quick Temp	✓	✓	✓

M3A

Vital Signs Monitor



Outpatient



Transport



Ward



Portable Design



Covidien Quick Temp



EDAN

A world of potential