

Viatom Checkme O2 MAX™



*Smart wrist oximeter with 0.5 Hz memory function
compatible with **Visi-Download** and with a **two year warranty***

The **Viatom Checkme O2 MAX™** is the next generation pulse oximeter for overnight monitoring of SpO₂ and pulse rate. Compatible with Visi-Download, the pulse oximeter uses a rechargeable battery and a revolutionary patented 'ring' pulse oximeter, for the best patient experience we have yet encountered.

- Four night recording memory
- Easy-to-clean, reusable ring sensor with **6 month warranty**
- Rechargeable battery
- Compatible with industry-standard Visi-Download (purchase may be required)
- Soft, comfortable silicone wrist strap, removable for cleaning
- Nail varnish? No problem! Sensor fits around base of thumb away from nails



The **Viatom Checkme O2 MAX™** is the next generation wrist-worn pulse oximeter.

Building upon the innovative design of the incredibly popular Checkme O2, Stowood have partnered with the manufacturer to bring **Visi-Download compatibility** to this brand new monitor.

With a simple wrist strap, single button, and combined download/charging cable, this device is easy to use and easy to clean.

The reusable ring sensor is the device's secret weapon - incredibly comfortable, you will probably forget you are wearing a monitoring device at all.





Small, light, with 72 hours' battery and room for four recordings.
Compatible with Visi-Download via the download/charging cable (included).

Patented ring sensor warranted for six months never comes off the patient's thumb, is comfortable and soft, and should lead to a better night's sleep than previous sensors (included).

This new Checkme O2 MAX uses the same sensors as the previous Checkme O2, and is backwards compatible with the Checkme O2 download cable and Visi-Download Viatom software download module (supplied separately). Standard Checkme O2 oximeters will still need the old download cable.

Technical specification

Recorded channels	SpO ₂ , pulse rate, motion
Dimensions	50 x 39 x 13 mm
Weight including batteries	25 g (main unit)
Power	Internal rechargeable batteries. 72 hours expected for typical recording parameters. Uses standard USB charger.
IP rating.....	IP22
Blood oxygen level range.....	70 – 99%
Pulse rate range	30 – 250 bpm
Blood oxygen accuracy	80 – 99%: ±2%, 70–79%: ±3%
Pulse rate accuracy.....	±2 bpm or ±2%
Sensors	Ring sensor included as standard. Soft tip sensor, paediatric soft sensor, single patient foam and neonatal non-woven SpO ₂ sensors available. Sensors are compatible with Checkme O2 and Checkme O2 MAX
Data storage	Internal memory holds four ten-hour recordings
Setup & signal check	Via OLED screen
Software	Visi-Download as used for oximeters and capnometers. Reports available in six languages
CE	Viatom holds CE 0197
Cleaning	Wrist strap and sensors are removable for cleaning via standard methods
Warranty	24 months from Stowood on main unit, 6 months on probes

viatom

Shenzhen Viatom Technology Co., Ltd.
 4E, Building 3, Tingwei Industrial Park, Baoan District, 518100 Shenzhen, P.R.China



© 2021 Stowood Scientific Instruments Ltd. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Stowood, their subsidiaries or affiliates.