

BT-200 Fetal Doppler

HI-bebe

HI-bebe is showing you what fetus is saying.

HI-bebe is pocket size doppler for Obstetricians, midwives and pregnant women.

- Compact and light
- High-Sensitivity Doppler probe
- Ergonomic design
- Large-Display & clean sound
- Easy to use
- Long time continuous use (6 hours)



Sound version



BT-200 Fetal Doppler **Hi-bebe**

Specifications

Display	• LCD Back light • Singnal quality display	• Low battery display
Heart rate range	• 50~240 bpm	
FHR Accuracy	• $\pm 2\%$	
Ultrasound Frequency	• 2 MHz	
Ultrasound Intensity	• $<10 \text{ m W/cm}^2$	
Speaker output power	• 1.2 W (Built in Loudspeaker)	
Audio Output for Earphone		
Auto Shut off Sound	• 1 Minute	
Auto Shut off Power	• 5 Minutes	
Power	• 1.5 V Battery X 2 (Type : AA)	
Sensitivity	• 10-12 Weeks onward	
PC Interface	• Sound card (BCM-220)	
Battery Life	• 360 Minutes	
Accessories	• Carrying case • Gel (60ml)	• Battery (1.5V AA X 2)

Waterproof (Optional)

PC Software

BCM220

- Record and play fetal heart sound (.wav format)

