

Electronic baby scale M112600

ADE

Large, but light-weight baby weighing scale ideal for GPs and for home visiting nurses as portable scale. The clear 30 mm display and the three key operations make the scale perfectly easy to use. Optionally available is a carrying bag with shoulder strap. Operated by standard AA batteries and mains.

Technical specification:

Capacity:	20 kg
Graduation:	5 g
Functions:	Tare, autohold, automatic zero setting, units switch kg/lbs, low battery indication
Dimensions:	600 x 280 x 40 mm
Weight approx.:	2,6 kg
Optional extra:	Carrying bag MZ10030

