



High-Definition **Multipurpose** Wireless Ultrasound

Image You Can Trust

Carry one scanner in your pocket with your smart phone or tablet to see high-resolution ultrasound imaging of deep structures and superficial anatomy.




Image Clarity You Can Trust for a Confident Diagnosis

Meet the world's only 3rd generation portable ultrasound. Now 30% smaller and lighter, Clarius HD3 delivers best-in-class portable ultrasound for your specialty, with an easy-to-use app powered by artificial intelligence and connected to the cloud.

- Ideal for abdominal, pelvic, and thoracic imaging, use Clarius for hip, cardiac, lung, OB/GYN, bladder, and superficial applications
- Scan within 20 seconds of starting the scanner and launching the Clarius App
- Choose from AI-powered presets for instant adjustments to get a high-definition image
- Tap to switch modes, pinch to zoom, slide up and down to change gain




Clarius C3 HD3 Ultra-High Frequency

Frequency: 2-6 MHz

Max Depth: 40 cm

Weight: 308g (0.68lb)


Dimensions: 146 x 76 x 32mm (5.7 x 3.0 x 1.3in)




Internal Jugular




Fetal Profile



Aortic-Aneurysm



Gallbladder



Lung



Subcostal



The images enable the ER team to quickly diagnose lung, heart, and abdominal concerns in our patients. During these difficult times, it has brought our team so much joy to receive such an amazing gift. We are so pleased with the new portable ultrasound and with how much more convenient and efficient it has already made our daily work."

Dr. Joanna Bostwick, MD - Emergency Physician

High-Resolution Ultrasound in a Portable Wireless Scanner



Clarius Cloud

Manage your exams and build reports anywhere, anytime with unlimited Clarius Cloud storage.



Easy To Use

Friendly AI-powered app for your Apple or Android automatically delivers the best imaging.



High-Definition

Octal Beamforming technology and 192 piezoelectric elements deliver high performance and resolution.



Whole Body Imaging

Versatile scanner provides superb imaging for the entire body.



Ultra-Portable

Ideal for small spaces, Clarius fits in your pocket. Take it anywhere!




Wireless

Steadfast connection. Total freedom. Without cords, cleaning and disinfection are simple and fast.




Since Clarius is so portable, I'm using ultrasound way more than I used to. It's like my stethoscope. The image quality is very good. With handheld ultrasound, we're now able to complete many ultrasound exams, including cardiac ultrasound, lung ultrasound, and abdominal ultrasound."

Dr. Virginia Robinson, MD - Family Physician




Clarius Membership


Clarius Membership empowers you with education options to hone your skills, advanced software tools that expand your capabilities, and a secure cloud solution that lets you take your practice anywhere. Get the Membership advantage to fast track your journey to mastering ultrasound.



Pulsed Wave Doppler



Auto Preset AI



Clarius Cloud

Additional Features

- Act One™
- Advanced Calcs
- Advanced Measurements
- Auto Heart Rate
- Clarius Live Telemedicine
- DICOM
- In-App Video Education
- Monthly Online Training
- OB Growth Tables
- Specialized Presets
- Split Screen Mode

Two Ways To Get Everything You Need

Membership Bundle

With your membership plan, you'll receive a discount on a scanner and the following:

- ✓ Support
- ✓ 60 Minutes of 1:1 Training
- ✓ Unlimited Clarius Cloud Storage
- ✓ Clarius Live
- ✓ Clarius Education
- ✓ Advanced Packages

One-Time Purchase

You may also purchase your scanner outright with the Advanced Packages à la carte:

- ✓ Support
- ✓ 60 Minutes of 1:1 Training
- \$ Advanced Packages à la carte

[REQUEST A QUOTE](#)

Clarius is on a mission to make accurate, easy-to-use, affordable ultrasound available to all medical professionals in every specialty. Our ultra-portable wireless scanners deliver best-in-class imaging with an easy-to-use app for iOS and Android that is powered by artificial intelligence and connected to the cloud. More than two million high-definition scans have been performed using Clarius wireless scanners, now available in over 90 countries worldwide.

Clarius Mobile Health

2985 Virtual Way #130,
Vancouver, BC, V5M 4X7
Canada
hello@clarius.com
Phone 1-778-800-9975

MKTG-00211 Rev 1