



EQUIPMENT CART

All EC models are well designed and high quality; with 4 castors 100 mm dia., 2 castors of which are lockable. The height of each shelf & each drawer is adjustable to fit your needs.



- Mount profile: aluminium
- Drawer/shelf: metal w/ powder coated, height adjustable
- Push handle: steel tube w/ powder coated
- Drawer/shelf cover: ABS
- Base cover: ABS
- Drawer rail: auto return



- Pole height adjustable by stepping
- Height range: 310mm (12")
- w/ push handle and 4" plastic casters
- Dimension:
W413 x D356 x H900~1210mm
(w/ 4" casters)
W523 x D550 x H900~1210mm
(w/ 4" casters & bumpers)



AEC4714

SPECIFICATIONS

EC5311 / EC5330 / EC5302 / EC5314

Cart dimensions (mm)	W500 x D610 x H1350 (w/ 4" casters) (EC5311)
Aluminium profile height (mm)	1230 / 1170 (above the base)
Drawer dimensions (mm)	W455 x D360 x H140
Shelf dimensions (mm)	W455 x D360 x H30
Cart weight (kg)	35 (EC5311)
Drawer load limit (kg)	10
Shelf load limit (kg)	10
Caster dimensions	4" steel casters (100mm, 2 w/ brake) Anti-static casters available at additional charge.



EVC

- VESA connector, height adjustable
- Material: steel
- Weight capacity: ~8kg (17.6 lbs)
- Rotation: 270°
- Pivot: 180°
- VESA: 75 / 100mm

MF SERIES EC



EIP + ELC iPad holder

ENH + ELC Notebook holder



ELC

- Monitor long arm, height adjustable
- Material: steel
- Max weight: 10kg
- Max length: 25cm
- VESA: 75 / 100mm

EKB

- Key board & mouse holder, height adjustable
- Tilt: 180°
- Material: aluminum alloy
- Extension: 540 mm



EHL

- Push handle, metal with powder coating, height adjustable

EPH

- Probe holder

EBH

- Bottle holder

EBB

- Utility wire basket, height adjustable (20 x 20 x 15cm)

